
# Flow Rate Sensor Specifications

### Media Isolated Microfluidic Flow Rate Sensor

- Measures up to 80 uL/min
- Non-invasive measurement
- Simple UART Output
- 15 Seconds Response Time



#### **Electrical Characteristics**

| Parameter      | Min | Nominal | Max     | Unit |
|----------------|-----|---------|---------|------|
| Supply Voltage | 1.8 | 3.3     | 3.6     | V    |
| Supply Current | -   | 10      | 600 (1) | mA   |
| Response Time  | 14  | 15      | 16      | s    |

Note 1: Drains up to 600 mA while applying power to the heater wires

#### **Measurement Characteristics**

| Parameter                    | Min | Nominal | Max | Unit         |
|------------------------------|-----|---------|-----|--------------|
| H2O Flow Rate Full Scale     | -   | -       | 80  | uL/min       |
| H2O Accuracy Below 55 uL/min | -   | -       | 5   | % full scale |
| H2O Accuracy Below 80 uL/min | -   | -       | 18  | % full scale |

All test results taken with the following parameters

- Supply Voltage = 3.3V
- Ambient Temperature = 23 C
- Water = DDW (Distilled, Deionized Water)

## Operation & Storage Characteristics

| Parameter                    | Min    | Nominal | Max    | Unit |
|------------------------------|--------|---------|--------|------|
| Storage Temperature          | -10    | 23      | 60 (2) | O    |
| Operating Temperature        | 10 (3) | 23      | 50 (3) | С    |
| Storage & Operating Humidity | -      | -       | 95     | %RH  |
| Maximum Reagent temperature  | -      | -       | 30 (4) | С    |

Note 2: The Tygon tubing may have a lower storage temperature than this figure

Note 3: These figures are best guesses, they have not been tested

Note 4: The maximum temperature the reagent will be heated to during the measurement process

### **Mechanical Characteristics**

| Parameter                           | Min   | Nominal | Max      | Unit |
|-------------------------------------|-------|---------|----------|------|
| Measurement Tube Inner Diameter (5) | 0.184 | -       | 0.254    | mm   |
| H2O Pressure Drop @ full scale      | -     | -       | 10 (6)   | mbar |
| Total Internal Volume               | -     | -       | 2        | uL   |
| Max Pressure                        | -     | -       | 2000 (7) | mbar |

Note 5: Measurement tube consists of alternating sections of stainless steel 29G tubes and 0.01"ID, 1/32"OD Tygon tubes

Note 6: Theoretical, calculated with 80uL/min as flow rate

Note 7: Experimental result where a force equivalent to 2 bars of pressure were applied to the syringe connected to the flow sensor


## **Specification Charts**

### Distilled, Deionized Water (DDW)


Benchmark Result with Fluigent System (SLI-0430)



#### Accuracy (%of f.s.)



#### Accuracy (%of m.v.)




|   | variable         | mean    | min  | median  | max     | nmissing | eltype  |
|---|------------------|---------|------|---------|---------|----------|---------|
| 1 | :abs_err         | 4.5975  | 1.6  | 3.25    | 14.15   | 0        | Float32 |
| 2 | :measurement_err | 14.8148 | 5.22 | 9.08187 | 65.4    | 0        | Float32 |
| 3 | :fullscale_err   | 5.74687 | 2.0  | 4.0625  | 17.6875 | 0        | Float32 |
| 4 | :flow            | 42.5    | 5    | 42.5    | 80      | Θ        | Int32   |

abs\_err: uncertainty in uL/min


measurement\_err: accuracy in %of measured value fullscale\_err: accuracy in %of full scale


flow: Data taken between these set flow rates

## Saline Solution (10% Concentration)




|   | variable         | mean    | min     | median | max  | nmissing | eltype  |
|---|------------------|---------|---------|--------|------|----------|---------|
| 1 | :abs_err         | 3.62    | 1.48    | 3.35   | 7.64 | 0        | Float32 |
| 2 | :measurement_err | 9.18634 | 2.58571 | 9.5375 | 17.0 | Θ        | Float32 |
| 3 | :fullscale_err   | 4.525   | 1.85    | 4.1875 | 9.55 | Θ        | Float32 |
| 4 | :flow            | 45.0    | 10      | 45.0   | 80   | Θ        | Int32   |







## 1x Phosphate Buffered Saline




|   | variable         | mean    | min     | median  | max     | nmissing | eltype  |
|---|------------------|---------|---------|---------|---------|----------|---------|
| 1 | :abs_err         | 5.1625  | 1.44    | 3.895   | 11.85   | 0        | Float32 |
| 2 | :measurement_err | 11.6892 | 6.23334 | 11.9083 | 14.8125 | 0        | Float32 |
| 3 | :fullscale_err   | 6.45312 | 1.8     | 4.86875 | 14.8125 | 0        | Float32 |
| 4 | :flow            | 45.0    | 10      | 45.0    | 80      | Θ        | Int32   |







### 10x Phosphate Buffered Saline



|   | variable         | mean    | min    | median  | max    | nmissing | eltype  |
|---|------------------|---------|--------|---------|--------|----------|---------|
| 1 | :abs_err         | 5.22125 | 2.43   | 5.095   | 7.89   | 0        | Float32 |
| 2 | :measurement_err | 14.1046 | 8.075  | 10.9881 | 35.7   | 0        | Float32 |
| 3 | :fullscale_err   | 6.52656 | 3.0375 | 6.36875 | 9.8625 | 0        | Float32 |
| 4 | :flow            | 45.0    | 10     | 45.0    | 80     | 0        | Int32   |




## Isopropyl Alcohol




Sensor requires recalibration to measure IPA.

## Reagent Temperature Variation


Water @ 10C




|   | variable         | mean    | min     | median  | max     | nmissing | eltype  |
|---|------------------|---------|---------|---------|---------|----------|---------|
| 1 | :abs_err         | 6.6025  | 3.38    | 6.2     | 10.53   | 0        | Float32 |
| 2 | :measurement_err | 18.4999 | 9.83334 | 13.3455 | 35.85   | 0        | Float32 |
| 3 | :fullscale_err   | 8.25313 | 4.225   | 7.75    | 13.1625 | 0        | Float32 |
| 4 | :flow            | 45.0    | 10      | 45.0    | 80      | 0        | Int32   |




#### Water @ 40C



|   | variable         | mean    | min    | median  | max    | nmissing | eltype  |
|---|------------------|---------|--------|---------|--------|----------|---------|
| 1 | :abs_err         | 6.2125  | 1.41   | 6.89    | 10.96  | 0        | Float32 |
| 2 | :measurement_err | 14.701  | 9.3125 | 14.3733 | 21.025 | 0        | Float32 |
| 3 | :fullscale_err   | 7.76562 | 1.7625 | 8.6125  | 13.7   | 0        | Float32 |
| 4 | :flow            | 45.0    | 10     | 45.0    | 80     | 0        | Int32   |



## **Pinout**



In the picture above, the pins are numbered 1 to 7 from bottom to top and 8 to 11 from left to right.

| Pin Number | Pin Description |
|------------|-----------------|
| 1, 9       | VDD             |
| 2          | SWCLK           |
| 3, 8       | GND             |
| 4          | SWDIO           |
| 5          | NRST            |
| 6          | VCP_TX          |
| 7          | VCP_RX          |
| 10,11      | Do not connect  |

#### Calibration

Manual calibration is required. Here are the resources required for the calibration step.

- 1. STLINK-V3
- 2. Access to the github repo https://github.com/jackyruth/Low\_Cost\_Microfluidic\_Flowsensor.git
- 3. A system to drive a reference flow rate

Here are the steps to calibrate the sensor

- 1. Connect the STLINK–V3 to the sensor, see the product specification sheet for the sensor pinout diagram
- 2. Program 'stm32/data collection' firmware
- 3. Drive a constant flow rate through the sensor and store the data. Preferably, the collected data corresponds to flow rates between 5 uL/min and 80 uL/min in increments of 5 uL/min
- 4. Push the data through 'julia/flowrate\_detection.jl' which outputs three calibration constants A, B and C
- 5. Change the calibration constants in 'julia/functional\_firmware'

Program the 'julia/functional\_firmware' and the device should start returning proper flow rate measurements.

### **Mechanical Specs**

Dimensions of the sensor are 28.75x53.00 mm.

The PCB is 0.3mm thin with conductively filled vias to reduce thermal resistance.

Mounting holes are for M2 screws.

The wetted materials are stainless steel tubing cut from 29G needles and 0.01" ID 1/32" OD Tygon tubing.